# AQUAJET®

## LD-A7

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## AQUAJET AND ORAL CARE

With regards to dental care, most people tend to focus on «whiteness» while ignoring the importance of gum tissue health. Failure to clean the surfaces of the teeth below the gum line, particularly between the teeth, aside from halitosis, results in the gradual accumulation of tartar and plaque. This accumulation allows ever increasing bacterial population, stimulating increased blood flow as the body's defenses are mobilized.

The first sign is often bleeding gums in response to whatever friction is applied. As the condition progresses, the gum tissue tends to move away from the tooth allowing the bacteria to penetrate further following the gum tooth interface, finally reaching and attacking the supporting bone structure. As the bone is eroded by bacterial action, the gums tend to fall in toward the receding bone. The result is the appearance of tooth elongation as more of the gum support and support structures fall back. As the tooth loses its root anchor, it gradually becomes loose and eventually falls out. Occasionally, the bacterial invasion of support bone structure can lead to abscesses, cysts and creeping infection threatening the brain itself.

Consequently, it cannot be over stressed the importance of maintaining scrupulous oral hygiene. The use of the AQUAJET after every meal is an important first step in a meticulous program of oral hygiene both for busy executive and his family.

The AQUAJET is endorsed and recommended by professional dentists. It cleans teeth crevices, between the teeth and under the gum line very effectively. The nozzle can be swiveled thus adjusting water direction and can access areas that brushing and flossing may not reach. The powerful water pressure flushes away food residue and massages the gums. It improves blood circulation and helps prevent periodontal disease, gingivitis, tooth decay and assists in inhibiting plaque formation. Using the device daily will promote and safeguard healthy teeth and gums. It is a remarkable, convenient tool for maintaining oral hygiene.

## IMPORTANT SAFEGUARDS

When using electrical products, especially when children are present, the basic safety precautions that are listed below should be always followed. Please read all instructions before using the AQUAJET.

### DANGER

- Always unplug this appliance immediately after using.
- Do not handle the plug with wet hands.
- Do not place or store product where it can fall or be pulled into a tub or sink, shower stall.
- Do not use while bathing.
- Do not place in or drop into water or other liquid.
- Do not reach for any appliance that has fallen into the water. Unplug it immediately.
- Before using, check the power cord for damage.

## WARNING

To reduce the risk of burns, electric shock, fire, serious injury or damage to persons:

- Do not use this appliance if the power cord or plug is damaged, if it is not working properly, if it has been dropped or damaged, or dropped into any liquid. Return the appliance to an authorized service center nearby for examination and repair. Do not try to fix it yourself.
- When the appliance is used by, on, or near children, close supervision is very necessary.
- An appliance should never be left unattended when it is plugged in. Use the appliance only for its intended use as described in this manual. Do not use attachments or accessories not recommended by the manufacturer.
- Keep the power cord away from hot or heated surfaces.
- Remove the Water Tank before filling. Fill with lukewarm water or appropriate oral rinse products only. Do not overfill.
- Never drop or insert any object into any opening or hose.
- Do not operate where oxygen or aerosol sprays are being used.
- Do not use outdoors.

## HOW TO USE

### Placing Water Tank

Take off the Water Tank or open the Cover of Water Tank to fill water or mouthwash water, if desired. Then, fit it back to its position (Fig.1). Make sure the water temperature does not exceed 40°C. Do not fill any hot water into the tank. The gum is very sensitive. It may be hurt by hot water.

### Power Cord And Plug

The plug and power cord is stored at the rear storage box (Fig. 2). Take off the rear panel by pushing it downward.

### Attaching/detaching nozzle to handle of irrigator

Attach one of the nozzles (choose the color you like) to the handle until it locks in place (Fig. 3).

#### **Recommendations for safe disconnection of nozzle:**

1. If the nozzle is detached from the handle badly or slide in the hand because of the water, use a dry cloth. Wrap the nozzle is not rough, dry cloth and repeat detaching.
2. When disconnecting the nozzle from the handle of irrigator drag it WITHOUT rotation around its axis. Rotation of the nozzle may lead to damage.

Please buy additional Nozzle LD-SA01 only.

### Power On Button

Press the power on button. For the first time use, please turn on the power on button for seconds until water coming out (Fig. 4).

*The use of the AQUAJET for dental hygiene cannot replace a toothbrush completely. Some friction is necessary. Use the AQUAJET after meals to flush away food particles and residues from the teeth and to provide a deep massage to the gums, following which use the toothbrush and dental floss as well to maximize oral hygiene and remove the plaque. Remember to floss before bedtime and irrigate vigorously afterwards with the AQUAJET; your teeth will love you for it!*

## Adjust Water Pressures

Adjust water pressure regulating knob before using (Fig 5). Turn it to the lowest setting (position 1) for first time use, especially if your gums tend to bleed or if children are using the AQUAJET. Gradually increase the water pressure over time until finding the most comfortable setting for yourself. Consult your professional dentists for recommendation of water pressure setting, particularly if the consumers have mouth disease.

## Start To Use

Bend over the nozzle tip and guide it to your teeth and gums. Slide the handle switch up to start the water stream (Fig. 6). Keep your mouth slightly open to allow the water to run into the washbasin. You can temporarily stop the water flow at any time by sliding down the switch of handle.

Direct the jet stream along the gum line with the jet tip. It can flush away the food particles remains on the teeth or between the teeth effectively and massage the gum gently.

## AFTER USE

According to studies of periodontal disease cases, the early stages of plaque formation are often ignored because of the lack of clear symptoms. Typically the disease is allowed to progress until the more serious middle or last stages. Bleeding from the gums upon use of the AQUAJET is symptomatic of periodontal disease, however, continued use of the AQUAJET tends to alleviate these symptoms and allows the gums return to a healthier state. If the bleeding carries on for more than two weeks, please consult your dentist.

After use, always empty the water tank, spray tube and nozzle entirely to prevent creating any bacteria from the standing water. Turn the power on button off. Take off the nozzle and store it into the nozzle retaining box. Wipe the handle dry and place it back onto its base.

## CLEAN AND MAINTENANCE

Before cleaning and maintaining the AQUAJET, always unplug from the electrical outlet first. Regularly clean the housing with a damp cloth. The Water Tank and Nozzles are removable and dishwasher-safe. When cleaning, please do not use any rough tool or utensil to wash the main body (such as scrubbing pad) to prevent the AQUAJET from damaging.

## CHECK AND INSTALLATION

Every unit is packed with 4 nozzles and a screw bag. After unpacking, please check them first. If there is any parts missing or broken, please contact your local dealer for replacement. The spray tube, nozzles or water tank may be with several drops of water left inside. It is left from the test. It is distilled water that is safe and hygienic. It is no harm to the health.

Do not repair or modify this AQUAJET yourself. If parts need to be replaced, changed or repaired, please contact your local dealer or an authorized service center.

## WALL MOUNTING INSTALLATION

Each unit of AQUAJET has been attached with one screw bag. The installed steps are listed as follows:

1. Take off the rear panel by pushing it downward.
2. Take out the screws from the bag.
3. Drill two holes of 0.8 and 2.5 deep. The distance between two holes is 8.3 cm.
4. Hammer down two plastic wall bolts and keep it in level with the wall.
5. Screw in the rear panel (Fig. 7).
6. Take the power cord out of the storage box and fit it through the hole of rear panel. Do not press too hard to damage the power cord (Fig. 8).
7. Please make sure to screw in tightly (Fig. 9).

## GUARANTEE

We grant a 1 year guarantee on the product commencing on the date purchase. The guarantee only comes into force if the date of purchase is confirmed by the dealer's stamp and signature on Guarantee Card.

## SPECIFICATIONS

S/N (LOT)	 LDA7
Model	LD-A7
Power source	200-240B AC ~50Hz
Power consumption	40W
Water pressure	29N~81N
Pump frequency	1200 imp./min.
Max. operating time	10 minutes
Water tank using time	1.5 minutes
Complete Set	Main Body, Water Tank, Rear Panel, 4 Nozzles, Screw Bag
Water Tank Volume	500 ml
External dimension	180(H) x 175(W) x 145(D)mm
Operating conditions:	
Temperature	+10°C to +35°C
Relative humidity	80% and below
Atmospheric pressure	86 to 106 kPa
Storage and transportation environment:	
Temperature	-10°C to +40°C
Relative humidity	95%
Atmospheric pressure	50 to 106 kPa
Net Weight	1350 g

## TROUBLE-SHOOTING GUIDE

Problem	Cause	Corrective Action
AQUAJET is no function.	«Power» is not connected.	Unplug it and plug it again.
	«Socket» has no electricity.	Try other sockets and plug it again.
Pump is functioning, but there is no water sprayed out	There is no water inside the water tank after using.	Refill water into the water tank. After seconds, the water stream will be jetted out.
Water pressure is too low.	It may be created due to the wrong water pressure regulation switch position.	Adjust the switch to the most suitable one.
	Water tank is placed wrong location	Check the water tank to see if it is placed to its correct location.

## ADDITIONAL ACCESSORIES\*

## Nozzle AQUAJET LD-SA01

- Intended for use with AQUAJET oral dental irrigator
- Made of plastic
- For individual use
- Numbers in packing – 2 pcs.

## Nozzle AQUAJET LD-SA02

- Subgingival nozzle for delivering therapeutic rinses into gum pockets, cleaning bridges and dentures.
- Intended for use with AQUAJET oral dental irrigator.
- Made of plastic with flexible nozzle.
- For individual use.
- Numbers in packing – 2 pcs.

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**LITTLE DOCTOR INTERNATIONAL (S) PTE. LTD.**